Initial Setup

You can open the app directly after installation or by tapping the app icon on your device. Then follow the Mobile Voice setup wizard.

Mobile Voice Setup Wizard

Complete the setup with your device in portrait orientation.

- 1. Open Mobile Voice by tapping the **app icon** on the home screen.
- Before you start going through the setup wizard a message asks if you want to allow Mobile Voice to access the microphone on your device. You need to allow this in order to record dictation. Tap Allow.



3. On Step 1 of the setup wizard, *read the instructions*, then tap the **Open Settings** link.



This takes you directly to the settings for Mobile Voice in your device Settings. Under ALLOW MOBILE VOICE TO ACCESS, tap Keyboards.

Mobile voice		< Apps	Mobile Voice	
Settings		ALLOW MOBI	ILE VOICE TO ACCESS	
Q. Search	Ų	U Micro	phone	
		Siri Siri		>
		Q Searc	ch	>
		🔲 Keybo	oards	>
· · · · ·		-		



5. Slide the toggle to enable **Keyboard**.



6. Once you enable that setting, the Allow Full Access setting appears. Slide the toggle to Allow Full Access



- 7. A dialogue pops up on screen asking for confirmation to allow full access for Mobile Voice. This includes a warning that full access allows nVoq to transmit anything you type, including sensitive information. *Mobile Voice is HIPAAcompliant, so that is OK.* Tap **Allow**.
- 8. Navigate back to Mobile Voice. Tap the **◄ Mobile Voice** link at the top-left of the screen.



< Mobile Voice	Mobile Voice Keyboards
Settings	Keyboard
Q Search	Allow Full Access
and the second second	When using one of these keyboards, the keyboard can access all the data you type. About Third-Party Keyboards & Privacy

9. Back in the setup wizard, tap the Next button.



Copyright 2000-2025, nVoq Incorporated. All Rights Reserved. | 20.0.0 V.1 06/13/2025



10. This advances the screen to Step 2: Switch to our Keyboard. *If your device is not already displaying the Mobile Voice keyboard*, **press and hold the Globe icon** at the bottom-left of the screen.



11. On the menu that appears, select Keyboard - Mobile Voice.

Keyboard Settings	4 6	5 6 7	8	9 0	
English (US)		Y U		<u>о</u> Р	delete
Emoji	D F	GH	JK	L	return
Keyboard — Mobile Voice _{English}	+ = V	/ ; B N	: ! M ;	?	shift
.?123				.?123	Ű

12. Your device should now be using the Mobile Voice keyboard. Tap the **GOT IT** button on the screen that displays.





- 13. You are taken to the sign in screen but must first respond to a popup message asking if you will allow Mobile Voice to track your activity.
 - Tap **Allow** if you want to be able to send crash reports to nVoq so we can fix issues in the keyboard that may occur.
 - Tap Ask App Not to Track if you do not want to send crash reports to nVoq.



Sign In Select Your System

On the sign in screen, leave the default system selected [Healthcare System (USA)].

TITLE BUILDED		
	Mobile Voice	
	Let's start transforming your clinical documentation!	
1. Select your system		
Healthcare System (USA)		~
2. Select your sign-in		
Continue with Microsoft S	SO	
O Continue with Okta SSO		
Sign in with nVoq usernam	le	
Having trouble deciding o	n a sign-in method? Contact your company administrator for assis	tance.

If you are using an iPhone 11 the sign in screen may look a little different, but the options are the same.



Sign In with nVoq Username

1. Tap the **Sign in with nVoq Username** button to sign in with that method.

TOTAL REPORT		
	Mobile Voice	
لا	Let's start transforming our clinical documentation!	
1. Select your system		
Healthcare System (USA)		~
2. Select your sign-in		
Continue with Microsoft SSO		
O Continue with Okta SSO		
Sign in with nVoq username		

2. When you tap this button you are taken to the nVoq sign in screen. Enter your **nVoq** username and **nVoq password**.

K Back		
	Mobile Voice	
	Let's start transforming your clinical documentation!	
Bob.burger@global.		
Password		⊳ ⊗
	Sign In	
		Forgot username or password?

Tap the eye- icon on the right side of the password field to see the password characters you type.



3. Tap the **Sign In** button.

	Mobile Voice	
	Let's start transforming your clinical documentation!	
Bob.burger@global.com		
Password.terrible.password		••••••••
	Sign In	
		Forgot username or password

- IF YOU ARE REQUIRED TO CHANGE YOUR PASSWORD ON INITIAL SIGN IN you will be taken directly to the Change Password where you can do that.. *Tap the eye icon to see the characters you type.*
- Once signed in the Mobile Voice dictation screen will be displayed. Swipe up from Home Indicator (the dark line at the bottom of the screen) to go back to the home screen where you can *select an application* into which you would like to dictate.





Use Dictation

- 1. Go to the application where you want to add your note.
- 2. Place your cursor in the text area where you want to add your note.
- 3. Tap the microphone / dictation button.



4. The dictation screen opens in an <u>active recording state</u>, which is indicated by the microphone icon at the bottom of the screen being red and no keyboard being present.









5. Say your dictation.



6. **Tap the red microphone icon** to stop recording. *Tapping anywhere on the screen will also stop your dictation.*





The microphone icon is black and the keyboard is present when Mobile Voice is no longer recording.



7. Optionally you can add to your dictation by *dictating (or typing)* on this screen. To dictate, tap the microphone icon again. (When this icon is red, Mobile Voice is recording.)

муемк		
	🧠 Mobile Voice	
This is an initial visit for diagnosis of CHF. He a long history of cardi fraction to approximat	or a 95-year-old gentleman with is referred to us by doctor Ande omyopathy and a recent decrea tely 20-25%.	a hospice rson, who reveals se in his ejection
The patient was seen clinical diagnosis of b that he has had signif 2-3 weeks. He is no le decreased and he had notably more weak. H living and he is increa night.	in the office approximately 3 da ronchitis. The FES, as well as th icant clinical and functional dec onger ambulatory. His appetite is d significant unmeasured weigh e is unable to perform any of the singly lowercase somnolent with	ys ago with a e daughter report line over the last s significantly t loss. He is e activities of daily n poor sleep at
All caps on physical e somnolent gentleman pressure is 90/60, he	xamination all caps off: Exam re n, who appears to be comfortable art rate is 80 and irregular	veals a cachectic e. His blood
pressure is 90/60, he	art rate is 80 and irregular	

You can place your cursor anywhere you want to add text and dictate (or type).



 Note: If you use the return key after a sentence that has no ending punctuation, the first word of the next sentence will not be capitalized automatically.



8. The send icon turns blue when the transcript is ready. **Tap the blue send icon** to send the transcript to your device clipboard.



9. After your tap the send icon, the transcript disappears from the dictation screen and is added to your clipboard. <u>You must go back to your application so that transcript will paste there</u>. If you immediately



start a new dictation, <u>the previous transcript will</u> <u>be lost</u>.

You can go back to your application by **tapping** the ◀ link at the top-left of the screen. (Alternatively you can swipe up from the black line at the bottom of the screen to view your open apps and then tap on your original application from there.)



10. Your transcript appears in your application where you last placed your cursor.

Once your transcript has been pasted into your application you can still use dictation to make changes:

- a. Select the text you want to change.
- b. Tap the dictation button
- c. On the dictation screen, say your new dictation.
- d. End your dictation, and then tap the Send icon.
- e. The selected text will be replaced by your new dictation.

You can also copy text from your application and edit it on the Mobile Voice dictations screen:

- a. Select the text you want to edit, then select "Copy" from the menu that appears.
- b. Tap the dictation button.
- c. Stop the dictation automatically in progress.
- d. Tap onto the dictation screen and select "Paste" from the menu that appears.
- e. The text appears on the dictation screen where you can make changes using dictation or typing.
- f. When you send the updated text back to your application it will include your changes and replace the text that you originally selected.

Tips for Better Dictation Accuracy

- Keep background noise to a minimum, if possible. When you're dictating, *audio quality* can make the difference between a good experience and a bad one. Many factors can influence audio quality, including ambient noise in the environment in which you're dictating.
- Speak at a normal, conversational rate. For best accuracy, be consistent in the pace and rhythm of your speech. It may be helpful to take a few seconds to think about what you want to say before you begin engage the microphone.
- If an acronym does not return correctly in a transcript, you can dictate the word "letter" before each letter to make individual letters return in dictations (for example, "*Letter* A *Letter* D *Letter* L"). If there are acronyms that you want to dictate frequently that do not return properly in your dictation transcript, they can be added to your dictation vocabulary by your nVoq account administrator.
- Use generic words (doctor, patient, etc.) instead of proper names whenever possible. Uncommon names are unlikely to be in the topic dictionary, and therefore unlikely to return successfully in the transcript. If there is a proper name that you want to dictate frequently and it does not return properly in your dictation transcript, it can be added to your dictation vocabulary. This can be done in nVoq Administrator by your nVoq account administrator.

While nVoq makes every effort to protect PHI, we recommend that users of any speech recognition technology follow the <u>HIPAA minimum necessary rule</u>, which states, "Protected health information (PHI) should not be used or disclosed when it is not necessary to satisfy a particular purpose or carry out a function".

<u>Not</u> dictating PHI will not only meet these federal requirements under HIPAA, but will also reduce concerns regarding poor accuracy that result from unique name spellings. In most cases if you are dictating into an EHR you have already identified the patient and do not need to dictate PHI.

ferrec	i to us	s by	doc	tor A	Ande	erso	n, w	hor	evea	ils a	long	his	tory	of ca	ardi	omy	opat	hy a	and a	recent
he pat ronchi inctior gnifica eak. H	tient v itis. TI nal de antly i le is u ise so	vas s he Fl cline decr unab	ES, a e ovi easi le to	in t as w er th ed a por ; with	he c rell a re la: nd h form h po	office is th st 2- ne ha n an ior s	e ap e da -3 w ad si y of leep	prov ugh eek gnif the at r	tern s. He ican actin	tely : epor e is n t unr vities	3 da t th no lo mea s of	nys a at hi nge sure daily	igo v e ha: r am ed wi y livii	rith a s had bula sight ng ar	a cli d sig itory t los nd h	nica mific t His s. H ie is	diag ant app incri	gno: clini petit hota basi	sis of ical a e is ibly n ngly	nd hore
EDIC/ e is or g eve tro pa latera	ATION n Core ry oth tch 12 I shou	IS: H eg 3. her d 2 hoi ulder	lis m 125 ay, L urs (pai	mg Lasio on ai n, ty	catic a da k 80 nd 1 pica	ns i ay, L mg 2 hc ally 3	inclu ynn twic ours 3-4 d	ide / so E ce di off. dose	Avelo Irill 5 aily, 1 He to Its pe	ox 4(i mg Inspr akes er da	00 n a di ra tv Tyl y,	ng, h ay, C vice enol	he is Coun daily I Extr	on d hadir , an a St	lay 4 n 2.5 d he reng	4 of 5 mg e is c gth e	a se a d in a f in a f	ven- ay, c tran (4-i	-day digoxi sderr 6 hou	dose. n 0.125 nal irs for
HYSIC	AL E	XAM	INA	TION	N: Ex	cam	reve	als	a ca	chec	tic	som	nole	nt ge	entle	ma	n, wi	no a	ppea	rs to be
omfort	Carries -	-													22-0-10				1.000	
omfort	table.	rii 5	0.00																	
mfort	labie.	riis																		
omfort	Labie.	riis																		
omfori	table.	rus																		
amfort	table.	rus																		
omfort	table.	rnia																		
omfort	a a	nis I	w		8	ſ	r				V				1	Ĩ	0		D	
ab	q		w	I	е	I	r	I	t	I	У	I	u	I	1	Ι	0	I	p	Gelete
ab caps loc	q	a	w	s	e	d	r	f	t	9	у	h	u	l	I	k	0	I	p	colote
tab caps loc	q x	a	w	s	e	d	r V	f	t	g	y	h	u M	I	1	k	0		p	costote return shift

Use the Keyboard

Next Word Predictive Typing Suggestions

The keyboard adapts to your typing by learning from previously entered text and displaying likely next-word suggestions in the bar above the keyboard.

1. Type a word on the keyboard, and then look to the bar above the keyboard to see suggestions for what

Patient	t																	
rep	orts			is				trie	ed			(can				was	
tab	q	W		е		r		t	I	у		u	I	i	I	0	р	delet
caps lo	_{ck} a		s	Ι	d	Ι	f		g	I	h	Ι	j		k	Ι	L	retur
shift	z	х	Ι	С	I	v		b	Ι	n		m	Ι	! ,	Ι	?		shit
	.?123						n)	/oq N	۸obil	e Voi	се						.?12	23

2. Tap on a word in the suggestion bar.



3. The word appears in your text.

Patien	was																		
		se	en										hos	pita	lizec	ł			Q
tab	q	w	Ι	е		r	I	t		у		u		i		0		р	delete
caps loo	_{ck} a		s	I	d	I	f	I	g		h		j		k		I		return
shift	z	х	Ι	с		v		b		n	Ι	m	I	! ,		?			shift
	.?123						n۷	oq N	lopil	e Voi	ce							.?123	Ĵ

- Next word suggestions are saved in a file as an account property every hour. A new file is downloaded once every 24 hours. These words are not saved automatically when you uninstall or upgrade the keyboard, but if you upgrade (or uninstall and reinstall), the last file that was saved is retrieved on sign in.
- Before words are added to the file they are checked against nVoq's suggestions database and Apple's words database. This is to ensure that words are valid so typos do not return as next word suggestions.

word may come next.

- The most recently used words that you typed after the previous word appear first in the suggestion bar, from newest to oldest.
- The list is not scrollable.
- On iPad, up to 5 suggestions display.
- On iPhone, up to 3 suggestions display.

Keyboard Functionality

You may see icons at the top of the keyboard (Undo, Redo, and Copy). *These are part of iOS, not Mobile Voice.* If you would like to remove them, go to **Settings > General > Keyboard > Deactivate "Shortcuts"**.



Switch between the QWERTY keyboard and the numbers/symbols keyboard, tap the 123 or ABC key.

			V
tab q w e r	t y u i o p	deleta tab 1 2 3 4 5 6 7	8 9 0 _{delete}
caps lock a s d	f g h j k l	return caps lock @ # \$ & * ()	, " return
shift Z X C V	b n m ! ?	shift #+= % - + = / ; :	! ? , . #+=
.2123	nVoq Mobile Voice .?12	23 🔲 🗰 ABC	ABC

When on the numbers and symbols keyboard you can **tap the #+= key** to use additional symbols. **Tap the 123 key** to return back to the regular numbers/symbols.

	Q														Q											
tab	1	2	3	4	5	6	7	8	9	0	delete	tab	1	2	3	4	5	6	7	8	9	0	delete			
caps lock	@		#	\$	&	·	()	•	"	return	caps lo	sk 🛃	E :	£	¥	-	^	[1	{	}	return			
#+=	%	-	+	=	1	;	:	1 7	?		#+=	.2123	§		~		1	<	>	1	?		.?123			
<u>م</u>	АВС	nVoq Mobile Voice ABC											ABC	nVoq Mobile Voice								ABC	Ũ			

To *capitalize a single letter* tap the **shift key** so that it turns white.



To *turn on caps lock* **tap the caps lock key** (or double-tap the shift key) so that the caps lock key turns white. To *turn caps lock off*, tap the **caps lock key** (or the shift key).



Double-tap the SPACE BAR to insert a period and a space (.) after the last character.

